


# Bottles for Injectables PE

20ml, 50ml, 100ml

**PLASTIC**  
Compatible with the environment



Ref. no	1679	1719	1680
Volume (ml)	20	50	100
Material	PE-HD	PE-HD	PE-HD
Colour	Natural	Natural	Natural
Dimension A (mm)	53.5	66.5	82.1
Dimension B (mm)	30.3	35.0	45.0
Dimension C (mm)	Ø32	Ø43	Ø52
Neck opening (mm)	Ø12.7	Ø12.7	Ø12.7
Weight (gr)	5.0	7.5	11.0
Brimful Volume (ml)	26.0	57.9	109.0
Packing (pcs/brick)	97	55	36
Packing (pcs/carton)	1,358	605	288
Volume of carton (m3)	0.094	0.094	0.094
Carton dimension l x w x h (mm)	600x400x391	600x400x391	600x400x391
Carton weight filled (kg)	8.7	5.5	4.7
Axial strength (filled 20°C) - (kg)	40	30	45
Cartons/EUR pallet	20	20	20
Quantity/EUR pallet	27,160	12,100	5,760
Sterilization	Beta ray		
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6		

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.